

Rev03
Update: Jun,09,2025

DATASHEET

KRAS WT Peptide (VVVGAGGVGK)

Cat. No.: RP30235

Overview

Description	Antibody screening, T-cell assays, Immune monitoring, Antigen specific T-cell stimulation, Cellular immune response
Sequence	{Val}{Val}{Val}{Gly}{Ala}{Gly}{Gly}{Val}{Gly}{Lys}
Sequence Shortening	VVVGAGGVGK
Molecular Weight	842.00

Properties

Purity	>95%
Solubility	Soluble in Ultrapure water under 1mg/ml
Form	Lyophilized
Storage	Store at -20°C. Keep tightly closed. Store in a cool dry place.

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China