

AAV9  
Update: Sep,04,2023**DATASHEET**

# AAV9

Cat. No.: RP30243CN

## Overview

<b>Description</b>	Pool of 182 peptides derived from a peptide scan (15mers with 11 aa overlap) through Capsid protein (Swiss-Prot ID: Q6JC40) of Adeno-Associated Virus (AAV) - Type 9.
<b>Sequence Shortening</b>	MAADGYLPDWLEDNLSEGIREWWALKPGAPQPKANQQHQDNARGLVLPGYKYLPGPNGLDKGEV NAADAAALEHDKAYDQQLKAGDNPYLKYNHADADEFQERLKEDTSFGGNLGRAVFQAKKRLLLEPLGL VEEAAKTAPGKKRPVEQSPQEPDSSAGIGKSGAQPAKKRLNFGQTGDTEVDPDPQPIGEPPAAPSGV GSLTMASGGGAPVADNNEGADGVGSSSGNWHCDSQWLGDREVTTSTRTWALPTYNNHLYKQISNS TSGGSSNDNAYFGYSTPWGYFDNRFHCHFSRDPWQRLINNNWGFPRKRLNFKLFNIQVKEVTDNN GVKTIANNLTSTVQVFTDSYQLPYVLGSAHEGCLPPFPADVMIPQYGYLTLNDGSQAVGRSSFYCL EYFPSQMLRTGNNFQFSYEFENVPFHSSYAHSQSLDRLMNPLIDQYLYLSKTINGSGQNQQTLKFS VAGPSNMAVQGRNYIPGPSYRQQRVSTTVTQNNNSEFAWPGASSWALNGRNSLMNPGPAMASHK EGEDRFFPLSGSLIFGKQGTGRDNVDADKVMITNEEEIKTTNPVATESYGVATNHQSAQAQAQTG WVQNQGILPGMVWQDRDVYLQGPWAKIPHTDGNFHPSPLMGGFGMKHPPPQILIKNTPVPADPP TAFNKDKLNSFITQYSTGQVSVEIEWELQKENSKRWNPEIQYTSNYKSNVFAVNTGVYSEPRPI GTRYLTRNL

## Properties

<b>Purity</b>	purity for each peptide sequence is 90%
<b>Form</b>	Lyophilized
<b>Storage</b>	Store at -20°C. Keep tightly closed. Store in a cool dry place.

生产商: 南京金斯瑞生物科技有限公司 江苏省南京市江宁区科学园雍熙路28号

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China