

Version: 02  
 Update: 09/30/2022

## DATASHEET

### pUC57 plasmid DNA

**Cat. No.** SD1176

**Name:** pUC57 plasmid DNA

**Size:** 10 µg


**Description:** pUC57 is a common used plasmid cloning vector in *E. coli*. The vector length is 2,710 bp and is isolated from *E. coli*

strain DH5α by standard procedures.

**Storage Buffer and Concentration:** Supplied in lyophilized form. Store at -20°C.

**Storage:** Store at -20°C.

### Map



**For research and manufacturing use. Direct human use, including taking orally and injection are forbidden.**

生产商: 南京金斯瑞生物科技有限公司 江苏省南京市江宁区科学园雍熙路28号

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28 Yongxi Road, Jiangning District, Nanjing, Jiangsu, China