


Mass Spectrum



Sample Information

User : Gary
 Month-Day Processed : 02/06/20
 Time Processed : 06:06:53 PM
 Injection Volume : 0.3
 Sample Name : Calcitonin Gene Related Peptide (CGRP) (8-37), human
 Cat No. : RP11089
 Lot No. : E6133920B
 Theoretical MW : 3125.59
 Observed MW : 3125.6

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment : GK11010007
 Interface Bias : +4.5kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2ml/min
 B.conc : 50%H2O/50%MeOH