


Mass Spectrum



Sample Information
Month-Day Processed : 09/28/20
Time Processed : 22:44:42
Injection Volume : 0.3
Sample Name : mTRP-2(180-188)
Cat No. : RP20250
Lot No. : E8011020I
Theoretical MW : 1175.34
Observed MW : 1174.3

Interface :ESI
Nebulizing Gas Flow :1.5L/min
CDL Temp :250
Block Temp :200

Equipment : GK11010007
Interface Bias : -3.5 kV
Drying Gas Flow :5 L/min
T.Flow :0.2 ml/min
B.conc :50%H2O/50%MeOH