


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 09/18/17  
Time Processed : 16:08:37  
Injection Volume : 0.2  
Sample Name : LSKL, Inhibitor of Thrombospondin (TSP-1)  
Cat No. : RPI9990  
Lot No. : PE6231709  
Theoretical MW : 459.58  
Observed MW : 459.45

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50% H<sub>2</sub>O/50% MeOH