


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 05/08/20  
Time Processed : 10:22:09 AM  
Injection Volume : 0.3  
Sample Name : HIF-1 alpha (556-574)  
Cat No. : RP20479  
Lot No. : E5746620E  
Theoretical MW : 2254.60  
Observed MW : 2253.6

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : -3.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H2O/50%MeOH