


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 04/12/18
Time Processed : 16:20:33
Injection Volume : 0.2
Sample Name : Jagged-1 (188-204)
Sample ID : U0260DD090-1
Theoretical MW : 2107.40
Observed MW : 2107.30

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH