


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 10/02/19
Time Processed : 06:01:39 PM
Injection Volume : 0.3
Sample Name : Kemptide (Phospho-Ser5)
Cat No. : RP20111
Lot No. : E7031819J
Theoretical MW : 851.89
Observed MW : 851.8

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200