

Sample Name : Gastrin-1, human  
 Cat. No. : RP12740  
 Lot. No. : PE5741610

Pump A : 0.065% trifluoroacetic in 100% water (v/v)  
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)  
 Total Flow:1 ml/min  
 Wavelength:220 nm


Time	Unit	Command	Value	Comment
0.01	Pumps	Pump A B.Conc	5	
25.00	Pumps	Pump A B.Conc	65	
25.01	Pumps	Pump A B.Conc	95	
31.00	Pumps	Pump A B.Conc	95	
31.01	Pumps	Pump A B.Conc	5	
40.00	Pumps	Pump A B.Conc	5	
40.01	Controller	Stop		

<<Column Performance>>

<Detector A>

Column : Alltima™ C18 4.6 x 250 mm

Chromatogram



1 Det.A Ch1 / 220nm

Peak Table

Detector A Ch1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	8.454	105574	13778	1.793
2	10.861	16515	2601	0.280
3	11.100	5731273	501734	97.335
4	12.850	5969	760	0.101
5	16.549	4939	585	0.084
6	17.948	18506	3183	0.314
7	20.888	5407	805	0.092
Total		5888184	523448	100.000