


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 10/27/16
Time Processed : 10:53:34
Injection Volume : 0.2
Sample Name : Gastrin-1, human
Cat No. : RPI2740
Lot No. : PE5741610
Theoretical MW : 2098.21
Observed MW : 2100.00

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50% H₂O/50% MeOH