


# Mass Spectrum



## Sample Information

Acquired by : Gary  
 Month-Day Processed : 11/28/17  
 Time Processed : 11:28:52  
 Injection Volume : 0.2  
 Sample Name : Competence-StimulatingPeptide-12261  
 Sample ID : U3265CK220-1  
 Cat. No. : RP20427  
 Lot. No. : PE0781711  
 Observed MW : 2154.00

Interface : ESI  
 Nebulizing Gas Flow : 1.5L/min  
 CDL Temp : 250°C  
 Block Temp : 200°C

Interface Bias : +4.5 kV  
 Drying Gas Flow : 5 L/min  
 T.Flow : 0.2 ml/min  
 B.conc : 50%H2O/50%MeOH