


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 08/12/17
Time Processed : 7:56:24
Injection Volume : 0.4
Sample Name : Bacterial Sortase Substrate III, Abz/DNP
Cat No. : RP20431
Lot No. : PE7121708
Theoretical MW : 928.00
Observed MW : 927.95

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH