


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 07/05/18
Time Processed : 12:59:55
Injection Volume : 0.2
Sample Name : HEL (46-61)
Cat No. : RP11736
Lot No. : PE1871807
Theoretical MW : 1753.83
Observed MW : 1754.10

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH