


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 11/25/17
Time Processed : 10:16:28
Injection Volume : 0.2
Sample Name : HIV-1 TAT Protein Peptide
Cat. No. : RP10129
Lot. No. : PE5301711
Theoretical MW : 1559.90
Observed MW : 1560.20

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias
Drying Gas Flow
T.Flow
B.conc

: +4.5 kV
: 5 L/min
: 0.2 ml/min
: 50% H₂O/50% MeOH