


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 01/09/20
Time Processed : 19:40:37
Injection Volume : 0.3
Cat No. : RP20140
Lot No. : E5291420A
Theoretical MW : 3787.20
Observed MW : 3786.5

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH

β -Amyloid (1-34)