


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 03/02/19
Time Processed : 11:46:37 AM
Injection Volume : 0.2
Sample Name : Magainin 2
Cat No. : RP11232-0.5
Lot No. : PE8431902
Theoretical MW : 2466.90
Observed MW : 2466.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH