


Mass Spectrum



Sample Information
Month-Day Processed : 11/29/20
Time Processed : 08:48:48 AM
Injection Volume : 0.3
Sample Name : Activated Protein C (390 - 404), human
Cat No. : RP20313
Lot No. : E2682120K
Theoretical MW : 1900.15
Observed MW : 1899.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H₂O/50%MeOH