

Sample Name :C-Type Natriuretic Peptide (CNP) (1-22), human
 Cat No. :RP11110
 Lot No. :E3443821K
 Time Processed : 16:08:58
 Month-Day-Year Processed : 11/29/2021

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm

Equipment: GK11010017

mV



1 Detector A Channel 1 / 220nm

Peak Table

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	4.302	26221	4940	0.682
2	11.252	29813	5366	0.776
3	11.599	10106	404	0.263
4	12.451	37167	1273	0.967
5	13.345	3465	339	0.090
6	13.621	24695	3871	0.643
7	13.785	3652078	490157	95.031
8	14.021	54070	6811	1.407
9	15.732	1540	331	0.040
10	16.536	3877	458	0.101
Total		3843033	513949	100.000