


Mass Spectrum



Sample Information

Acquired by: Gary
Month-Day Processed: 06/20/19
Time Processed: 10:04:27 AM
Injection Volume: 0.4
Sample Name: LSKL, Inhibitor of Thrombospondin (TSP-1)
Cat. No.: RP19990
Lot. No.: E1187919F
Theoretical MW: 459.58
Observed MW: 459.8

Interface: ESI
Nebulizing Gas Flow: 1.5L/min
CDL Temp: 250
Block Temp: 200

Equipment: GK11010007
Interface Bias: +4.5 kV
Drying Gas Flow: 5 L/min
T.Flow: 0.2 ml/min
B.conc: 50% H₂O/50% MeOH