


Mass Spectrum



Sample Information

User :	Gary	Equipment	: GK11010007
Month-Day Processed :	02/15/20	Interface Bias	: +4.5kV
Time Processed :	18:53:15	Drying Gas Flow	: 5 L/min
Injection Volume :	0.3	T.Flow	: 0.2ml/min
Sample Name :	Angiotensin I, human	B.conc	: 50%H2O/50%MeOH
Cat No. :	RP11251		
Lot No. :	E6378220B		
Theoretical MW :	1296.48		
Observed MW :	1296.3		

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200