


# Mass Spectrum



## Sample Information

Acquired by : Gary  
Month-Day Processed : 07/22/20  
Time Processed : 20:09:30  
Injection Volume : 0.5  
Sample Name : Bombesin  
Cat. No. : RP11157  
Lot. No. : E9668220G  
Theoretical MW : 1619.90  
Observed MW : 1619.4

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% H2O/50% MeOH