


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 10/03/19  
Time Processed : 07:19:53 AM  
Injection Volume : 0.3  
Sample Name : HIV p17 Gag (77-85)  
Cat No. : RP17316  
Lot No. : E7025519J  
Theoretical MW : 981.11  
Observed MW : 981.0

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH