


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 09/25/19
Time Processed : 08:52:37 PM
Injection Volume : 0.3
Sample Name : PTD-p65-P1 Peptide
Cat No. : RP20522
Lot No. : E4942719I
Theoretical MW : 3829.50
Observed MW : 3828.8

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH