


Mass Spectrum



Sample Information

Acquired by : Gary
 Month-Day Processed : 07/25/18
 Time Processed : 15:43:47
 Injection Volume : 0.2
 Sample Name : Gastrin-1, human
 Cat. No. : RP12740
 Lot. No. : PE8301807
 Theoretical MW : 2098.21
 Observed MW : 2099.6

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250°C
 Block Temp : 200°C

Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50% H₂O/50% MeOH