


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 01/05/18
Time Processed : 12:28:58
Injection Volume : 0.4
Sample Name : LSKL, Inhibitor of Thrombospondin (TSP-1)
Cat No. : RPI9990
Lot No. : PE2721801
Theoretical MW : 459.58
Observed MW : 459.40

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH