


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 12/01/19
Time Processed : 23:40:56
Injection Volume : 0.4
Sample Name : Tetanus toxin (830-843)
Cat No.: RP19980
Lot No.: E7088819K
Theoretical MW : 1611.84
Observed MW : 1611.6

Interface : ESI

Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH

: GK11010007