


Mass Spectrum



Sample Information

Acquired by: Gary
 Month-Day Processed: 03/01/17
 Time Processed: 13:15:27
 Injection Volume: 0.2
 Sample Name: IRBP (1-20), human
 Cat No.: RP20269
 Lot No.: PE0341702
 Theoretical MW: 2194.60
 Observed MW: 2194.30

Interface: ESI
 Nebulizing Gas Flow: 1.5L/min
 CDL Temp: 250°C
 Block Temp: 200°C

Interface Bias: +4.5 kV
 Drying Gas Flow: 5 L/min
 T.Flow: 0.2 ml/min
 B.conc: 50%H2O/50%MeOH