


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 11/12/16
Time Processed : 21:30:45
Injection Volume : 0.2
Sample Name : GP (33-41)
Cat No. : RP20091
Lot No. : PE4831611
Theoretical MW : 1016.18
Observed MW : 1016.10

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH