


Mass Spectrum



Sample Information

User : Gary
 Month-Day Processed : 03/15/19
 Time Processed : 11:02:04 AM
 Injection Volume : 0.2
 Sample Name : C-Type Natriuretic Peptide (CNP) (1-22), human
 Cat No. : RP11110
 Lot No. : PE5681903
 Theoretical MW : 2197.60
 Observed MW : 2197.2

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment : GK11010007
 Interface Bias : +4.5kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2ml/min
 B.conc : 50%H2O/50%MeOH