


# Mass Spectrum



Sample Information  
Month-Day Processed : 12/29/20  
Time Processed : 8:18:42 PM  
Injection Volume : 0.3  
Sample Name : Angiotensin I, human  
Cat No. : RP11251  
Lot No. : E2662920L  
Theoretical MW : 1296.48  
Observed MW : 1296.2

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +3.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H2O/50%MeOH