


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 08/05/18
Time Processed : 18:23:48
Injection Volume : 0.2
Sample Name : Alpha 1(I) Collagen (614-639), human
Cat No. : RP20424
Lot No. : PE6071808
Theoretical MW : 2942.20
Observed MW : 2942.5

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250 °C
Block Temp : 200 °C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH