


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 03/03/20
Time Processed : 05:48:14 PM
Injection Volume : 0.6
Sample Name : Calcitonin Gene Related Peptide (CGRP), rat
Cat No. : RP11095
Lot No. : E9474920B
Theoretical MW : 3806.25
Observed MW : 3805.2

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200