


Mass Spectrum



Sample Information

User :	Gary	Equipment	: GK11010007
Month-Day Processed :	09/18/19	Interface Bias	: +4.5kV
Time Processed :	03:08:58 PM	Drying Gas Flow	: 5 L/min
Injection Volume :	0.3	T.Flow	: 0.2ml/min
Sample Name :	C-telopeptide	B.conc	: 50%H ₂ O/50%MeOH
Cat No. :	RP10566		
Lot No. :	E4746719I		
Theoretical MW :	868.90		
Observed MW :	868.8		

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200