


# Mass Spectrum



## Sample Information

Acquired by : Gary  
Month-Day Processed : 02/20/17  
Time Processed : 13:32:41  
Injection Volume : 0.2  
Sample Name : Proctolin  
Sample ID : U1633CB130-1  
Theoretical MW : 648.75  
Observed MW : 648.40

## Interface

:ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

## Interface Bias

: +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% H<sub>2</sub>O/50% MeOH