


Mass Spectrum



Sample Information
Month-Day Processed : 07/31/21
Time Processed : 9:10:03
Injection Volume : 0.2
Sample Name : Bombesin
Cat. No. : RP11157
Lot. No. : E8358021G
Theoretical MW : 1619.90
Observed MW : 1620.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : ZJ21010035
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH