


Mass Spectrum



Sample Information
Month-Day Processed : 04/19/21
Time Processed : 7:31:23
Injection Volume : 0.4
Sample Name : LCMV gp33 - 41
Cat No : RP20257
Lot No : E5236421D
Theoretical MW : 1044.25
Observed MW : 1043.5

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200
Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH