


Mass Spectrum



Sample Information
User : Gary
Month-Day Processed : 04/02/20
Time Processed : 02:42:58 PM
Injection Volume : 0.3
Sample Name : C-Peptide (3-33), human
Cat No. : RP11116
Lot No. : E8770420D
Theoretical MW : 3020.26
Observed MW : 3020.4

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200