


Mass Spectrum



Sample Information

Month-Day Processed : 03/20/21
Time Processed : 4:32:23 PM
Injection Volume : 0.5
Sample Name : HIF-1 alpha (556-574)
Cat. No. : RP20479
Lot. No. : E0991421C
Theoretical MW : 2254.60
Observed MW : 2253.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200
Equipment : GK11010007
Interface Bias : -3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH