

Sample Name :Peptide YY (PYY) (3-36), human  
 Cat. No. : RP10354  
 Lot. No. : PE8151902  
 Time Processed : 8:35:23  
 Month-Day-Year Processed : 03/06/2019

Pump A : 0.065% trifluoroacetic in 100% water (v/v)  
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)  
 Total Flow:1 ml/min  
 Wavelength:220 nm

Time	Unit	Command	Value	Comment
0.01	Pumps	Pump B Conc.	5	
25.00	Pumps	Pump B Conc.	65	
25.01	Pumps	Pump B Conc.	95	
32.00	Pumps	Pump B Conc.	95	
32.01	Pumps	Pump B Conc.	5	
40.00	Pumps	Pump B Conc.	5	
40.01	Controller	Stop		


<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm

Equipment: G09010878

Chromatogram



1 Det.A Ch1 / 220nm

Peak Table

Detector A Ch1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	13.283	8216	785	0.050
2	13.566	51955	7862	0.319
3	13.783	355723	64720	2.182
4	13.915	15522672	1548303	95.234
5	14.175	302517	41979	1.856
6	14.566	58477	4360	0.359
Total		16299560	1668009	100.000