


Mass Spectrum



Sample Information

Acquired by : Gary
 Month-Day Processed : 12/15/19
 Time Processed : 04:44:10 PM
 Injection Volume : 0.8
 Sample Name : TAT peptide
 Cat. No. : RP20256
 Lot. No. : E5858219L
 Theoretical MW : 1621.91
 Observed MW : 1622.0

Interface

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment

Equipment : GK11010007
 Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50% H2O/50% MeOH

Cat. No. : RP111133