


Mass Spectrum



Sample Information

Month-Day Processed : 05/15/21
Time Processed : 18:40:39
Injection Volume : 0.2
Sample Name : HEL (46-61)
Cat. No. : RP11736
Lot. No. : E2678521E
Theoretical MW : 1753.83
Observed MW : 1753.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH