

Sample Name :Peptide YY (PYY) (3-36), human

Cat. No. : RP10354

Lot. No. : E0334719F

Time Processed :5:57:14

Month-Day-Year Processed :06/14/2019

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

Time	Unit	Command	Value	Comment
0.01	Pumps	Pump B Conc.	5	
25.00	Pumps	Pump B Conc.	65	
25.01	Pumps	Pump B Conc.	95	
27.00	Pumps	Pump B Conc.	95	
27.01	Pumps	Pump B Conc.	5	
35.00	Pumps	Pump B Conc.	5	
35.01	Controller	Stop		


<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm

Equipment :SS-CM-0310

Chromatogram



1 Det.A Ch1 / 220nm

Peak Table

Detector A Ch1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	13.542	4199	789	0.064
2	15.000	9700	1008	0.148
3	15.132	11024	1508	0.168
4	15.258	10948	1683	0.167
5	15.602	58072	3525	0.886
6	16.189	6448718	646650	98.361
7	17.426	5536	1007	0.084
8	17.769	1405	207	0.021
9	18.144	4103	829	0.063
10	20.068	2493	295	0.038
Total		6556199	657500	100.000