


# Mass Spectrum



## Sample Information

Month-Day Processed : 11/13/20  
 Time Processed : 12:17:07 AM  
 Injection Volume : 0.5  
 Sample Name : Exendin-4  
 Cat. No. : RP10874  
 Lot. No. : E4703320K  
 Theoretical MW : 4186.66  
 Observed MW : 4185.6

## Interface : ESI

Nebulizing Gas Flow : 1.5L/min  
 CDL Temp : 250  
 Block Temp : 200

## Equipment

Interface Bias : +4.5 kV  
 Drying Gas Flow : 5 L/min  
 T.Flow : 0.2 ml/min  
 B.conc : 50% H<sub>2</sub>O/50% MeOH

Equipment : GK11010007