


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 06/21/18
Time Processed : 20:38:29
Injection Volume : 0.2
Sample Name : PACAP (1-27), human, ovine, rat
Cat No. : RPI0334
Lot No. : PE9941806
Theoretical MW : 3147.62
Observed MW : 3148.00

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH