

Sample Name :Glucagon-Like Peptide (GLP) I (7-37)  
 Cat No.: RP12738  
 Lot No.:E0603819L  
 Time Processed :8:58:01  
 Month-Day-Year Processed :01/09/2020

Pump A : 0.065% trifluoroacetic in 100% water (v/v)  
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)  
 Total Flow:1 ml/min  
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	25
25.00	Pumps	Pump A B.Conc	85
25.01	Pumps	Pump A B.Conc	95
32.00	Pumps	Pump A B.Conc	95
32.01	Pumps	Pump A B.Conc	25
40.00	Pumps	Pump A B.Conc	25
40.00	Controller	Stop	


<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm  
 Equipment:ZJ19010324

### <Chromatogram>

mV



### <Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	12.245	4657	518	0.114
2	14.418	7106	599	0.174
3	15.523	67083	6967	1.640
4	15.914	3912040	383314	95.618
5	16.389	73605	5519	1.799
6	17.283	10595	1748	0.259
7	17.933	6075	1004	0.148
8	18.347	3550	666	0.087
9	24.413	6597	802	0.161
Total		4091308	401138	100.000