


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 03/26/19
Time Processed : 09:11:26 PM
Injection Volume : 0.2
Sample Name : Brain Natriuretic Peptide (BNP) (1-32), human
Cat No. : RP11119
Lot No. : PE2251903
Theoretical MW : 3464.10
Observed MW : 3464.4

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH