


# Mass Spectrum



## Sample Information

Acquired by : Gary  
Month-Day Processed : 01/17/18  
Time Processed : 16:47:56  
Injection Volume : 0.2  
Sample Name : Endothelin 1, human, porcine  
Cat. No. : RP12758  
Lot. No. : PE1571801  
Theoretical MW : 2491.90  
Observed MW : 2492.10

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% H2O/50% MeOH