


# Mass Spectrum



## Sample Information

Month-Day Processed : 01/23/21

Time Processed : 3:37:24 PM

Injection Volume : 0.4

Sample Name : HIV-1 TAT Protein Peptide

Cat. No. : RP10129

Lot. No. : E7137521A

Theoretical MW : 1559.90

Observed MW : 1560.0

## Interface

Nebulizing Gas Flow : 1.5L/min

CDL Temp : 250

Block Temp : 200

## Equipment

Interface Bias : +4.5 kV

Drying Gas Flow : 5 L/min

T.Flow : 0.2 ml/min

B.conc : 50% H<sub>2</sub>O / 50% MeOH

Equipment : GK11010007