


# Mass Spectrum



Sample Information  
User : Gary  
Month-Day Processed : 09/19/19  
Time Processed : 7:30:05  
Injection Volume : 0.3  
Sample Name : LL37 (Human)  
Cat No. : RP13323  
Lot No. : E4880919I  
Theoretical MW : 4493.28  
Observed MW : 4493.4

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH